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MONTREAL
TRADING COMPANY

**Per-Oxy Activated
Hydrogen peroxide...0.5% (w/w)**

DIN 02504863

**BACTERICIDE
VIRUCIDE
FUNGICIDE**

**READY-TO-USE
HARD SURFACES CLEANER AND DISINFECTANT
RTU**

DESCRIPTION

Per-Oxy Activated is a disinfectant cleaner based on activated Hydrogen Peroxide technology to deliver fast and effective

cleaner performance.

- Highly Effective, ready-to-use disinfectant
- Disinfects hard surfaces as a virucide, bactericide and fungicide
- For use as cleaner and disinfectant of hard, nonporous surfaces

USE DIRECTIONS

For Hospitals, Healthcare and Dental clinics:

This product is not to be used as a sterilant/high-level disinfectant on any surface or instrument that: (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Disinfection of non-critical items, equipment and hard surfaces coming in contact with intact skin such as hemodialysis machines, stethoscopes, tabletops etc.: For heavily soiled surface, a pre-cleaning step is required; apply *as is* to surfaces with a sprayer, cloth or disposable wipe. Ensure surface remains wet for 10 minutes. Then, wipe surface dry or rinse.

For Industrial/institutional and Food processing: For heavily soiled surface, a pre-cleaning step is required; apply *as is* to surfaces and equipment; allow to remain wet for 10 minutes. Then, wipe dry. No rinsing is required.

For food-processing establishments: Food contact surfaces require rinsing with potable water. Avoid contamination of food during application and storage

For Barns:

-Remove all animals/poultry and their feed from premises prior to disinfection;

-Remove all heavy soil, such as urine and fecal matter, from hard surfaces and objects prior to disinfection;

-Empty all feeding and watering appliances prior to disinfection; and

-Following disinfection, do not house animals/poultry until areas have been ventilated.

For heavily soiled surface, a pre-cleaning step is required; apply *as is* to surfaces and equipment; allow to remain wet for 10 minutes. Then, wipe dry. No rinsing is required. **All**

surfaces and/or objects that will contact feed or drinking water should be rinsed with potable water before reuse.

PROPERTIES

Appearance: Colorless liquid / Odour: Odorless

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